

sequence to be used as a colony primer can include any sequence, but preferably includes 5'-AGAAGGAGAAGGAAAGGGAAAGGG (SEQ ID NO:1) or 5'-CACCAACCCAAACCAACCCAAACC (SEQ ID NO:2). The colony primer can be 5 to 100 bases in length, but preferably 15 to 25 bases in length. Naturally occurring or non-naturally occurring nucleotides may be present in the primer. One or two different colony primers may be used to generate nucleic acid colonies in the methods of the present invention. The colony primers for use in the present invention may also include degenerate primer sequences.--

Please replace the paragraph beginning at the bottom of page 44, line 28, with the following rewritten paragraph:

-- Figure 8: shows hybridization of probes to oligonucleotides attached to Nucleolink, before and after PCR cycling. The figure shows R58 hybridization to CP2 (5'-(phosphate)- TTTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 closed circles, CP8 (5' (amino-hexamethylene)- TTTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 closed triangles, CP9 (5' (hydroxyl)- TTTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 diamonds, CP10 (5' (dimethoxytrityl)-TTTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 open circles and CP11 (5' (biotin)- TTTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 open triangles.--

Please replace the paragraph beginning at page 45, line 16, with the following rewritten paragraph:

--The properties of the colony primers have been chosen based on a selection for oligonucleotide primers that show little non-specific nucleotide incorporation in the presence of heat-stable DNA polymerases. The colony primers, CP (5'-p CACCAACCCAAACCAACCCAAACC) SEQ ID NO:2 and CP (5'-p AGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:1 have been selected due to their low incorporation of radiolabeled [^{32}P -dCTP] in the presence of a stable DNA polymerase (AmpliTaq, Perkin Elmer, Foster City, CA) in the standard buffer and under thermocycling conditions (94 C for 30 seconds, 65 C for 1 minute, 72 C for 2 minutes, 50 cycles).--

Please replace Table 1, beginning on page 48, with the following new Table 1:

-- TABLE I

List of oligonucleotides used for templates preparation and colonies generation:

Name	DNA sequence	Coordinates (orientation)	Oligonucleotide Modification	Use
TP1	GAGGCCAGAACAGTTCAAGG (SEQ ID NO:3)	9810 (R)		Template 3.2 Kb
TP2	CCTGTGACAAGACGACTGAA (SEQ ID NO:4)	6550 (F)		Template 3.2 Kb
CP1	TTTTTTTTTTCACCAACCCAAAC CAACCCAAACC (SEQ ID NO:5)	None	5'P	Generate colonies
CP2	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:6)	None	5'P	Generate colonies
CP3	TTTTTTTTTTCACCAACCCAAAC CAACCCAAACC (SEQ ID NO:7)	None	5'SH	Generate colonies
CP4	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:8)	None	5'SH	Generate colonies
CP5	AGAAGGAGAAGGAAAGGGAAGG GTTTTTTTTTTTTTTNN (SEQ ID NO:9)	None	5'P	Generate colonies
CP6	AGAAGGAGAAGGAAAGGGAAGG GGG (SEQ ID NO:10)	None	5'P	Generate colonies
CP7	TTTTTTTTTTCACCAACCCAAAC CAACCCAAACC (SEQ ID NO:5)	None	5' (NH ₂)	Generate colonies
CP8	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:6)	None	5' (NH ₂)	Generate colonies
CP9	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:6)	None	5' (OH)	Control oligo
CP10	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:6)	None	5' (DMT)	Control oligo
CP11	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGG (SEQ ID NO:6)	None	5' (biotin)	Control oligo
TPA1	AGAAGGAGAAGGAAAGGGAAGG GCCTGTGACAAGACGACTGAA (SEQ ID NO:12)	6550 (F)	5'P	Template A
TPA2	TTTTTTTTTTAGAAGGAGAAGGA AAGGGAAGGGCGGCCGCTGAG GCCAGTGAAGTCAGA (SEQ ID NO:13)	7403 (R)	5'P	Template A
<u>TPB3</u>	TTTTTTTTTTCACCAACCCAAAC CAACCCAAACCGAGCTCAGGCTG AGGCAGGAGAATTG (SEQ ID NO:14)	9049 (F)	None	Template B'
TPB1	AGAAGGAGAAGGAAAGGGAAGG GGAGCTGAGGAGGAAGAGAGG (SEQ ID NO:15)	9265 (F)	None	Template B
TPB2	AGAAGGAGAAGGAAAGGGAAGG GGCGGCCGCTCGCCTGTTCTGG AAGACA (SEQ ID NO:16)	8411 (R)	5'P	Template B
TPB4	AGAAGGAGAAGGAAAGGGAAGG GGCGGCCGCTCGCCTGTTCTGG AAGACA (SEQ ID NO:11)	9265 (R)	5'SH	Template B'

Coordinate from HUMOXRAGE gene Accession number D28769
(R) means "reverse" and (F) means "forward"--

Please replace the paragraph beginning at the bottom of page 49, line 15, with the following rewritten paragraph:

-- A 3.2 Kb DNA fragment was taken as a model system to demonstrate the feasibility of colony generation from random primer PCR amplification. This strategy can be applied to sequencing of DNA fragments of approximately 100 Kb in length and, by combination of fragments to whole genomes. A fragment of DNA of 3.2 Kb was generated by PCR from human genomic DNA using PCR primers; TP1 5'-pGAGGCCAGAACAGTTCAAGG (SEQ ID NO:3) and TP2 5'-pCCTGTGACAAGACGACTGAA (SEQ ID NO:4) as described in example 1. The 3.2 Kb fragment was cut in smaller fragments by a combination of restriction enzymes (EcoRI and HhaI yielding 4 fragments of roughly 800 bp). The cut or uncut fragment DNAs were then mixed with the degenerate primer, p252 (5'-P TTTTTTTTTTISISISISISIS, SEQ ID NO:17 where I stands for inosine (which pairs with A, T and C) and S stands for G or C) and covalently coupled to the Nucleolink wells (Nunc, Denmark). The tubes were then subjected to random solid phase PCR amplification and visualized by hybridisation with labeled DNA probes, as will be described in Example 2a.--

Please replace the paragraph beginning at page 50, line 17, with the following rewritten paragraph:

-- A colony primer (CP2, 5'-TTTTTTTTTTAGAAGGAGAAGGAAAGGGAAAGGG) (SEQ ID NO:6, phosphorylated at its 5' terminus (Microsynth GmbH, Switzerland), was attached onto Nucleolink plastic microtitre

wells (Nunc, Denmark) in the presence of varying doses of Template A (prepared as described in example 1). 8 wells were set up in duplicate with seven 1/10 dilutions of template with CP2, starting with the highest concentration of 1 nM.--

Please replace the paragraph beginning at page 51, line 24, with the following rewritten paragraph:

-- Colonies Visualization

Probe: The probe was a DNA fragment of 1405 base pairs comprising the sequence of the template at their 3' end (nucleotide positions 8405 to 9259). The DNA probe was synthesized by PCR using two primers: p47 (5'-GGCTAGGAGCTGAGGAGGAA) SEQ ID NO:20, amplifying from base 8405, and TP2, biotinylated at 5' end, amplifying from base 9876 of the antisense strand.--

Please replace the paragraph beginning at page 54, line 21, with the following rewritten paragraph:

--A colony primer (CP2:

5'pTTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:8, phosphorylated at its 5' termini (Microsynth GmbH, Switzerland), was grafted onto Nucleolink plastic microtitre wells (Nunc, Denmark) in the presence of varying doses of the two templates A and B (prepared as described in example 1). Series of 8 wells were set up in triplicate with seven 1/10 dilutions of both templates starting with the highest concentration of 1 nM. Template dilutions are set up in opposite directions such that the highest concentration of one

template coincides with the lowest of the other.--

Please replace the paragraph beginning at page 67, line 5, with the following rewritten paragraph:

-- Colony primers CP1

(5'-pTTTTTTTTTTCACCAACCCAAACCAACCCAAACC) SEQ ID NO:7 and CP2
(5'-pTTTTTTTTTTAGGAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 which are
5'phosphorylated (Microsynth GmbH, Switzerland) and DNA
template B (prepared as described in example 1) were 5'
covalently attached onto 5 mm diameter glass slides (Verrerie
de Carouge, Switzerland) to a final concentrations of 1 μ M and
10 nM respectively, as follows: 2 nmoles of each primer were
added to 0.2 nmoles of template in 1 ml of solution A (41 μ l
of Methylimidazole (Sigma, #M-8878) in 50 ml H₂O, pH adjusted
to 7 with HCl) and then mixed 1:1 with solution D (0.2 mM EDC
in 10 ml of solution A). On both glass slides sides, 3.5 μ l
of the mixture were loaded, and incubated over night at room
temperature. The glass slides were then briefly rinsed with
5xSSC buffer and placed at 100 C in 10mM Tris buffer pH 8.0
for 2x5'.--

Please replace the paragraph beginning at page 67, line 27, with the following rewritten paragraph:

-- Glass slides were then individually placed onto a
MicroampTM reaction tube (Perkin Elmer) containing 170 μ l of
PCR mix, and DNA colonies were then generated using Taq
polymerase (AmpliTaq, PE-Applied Biosystems Inc., Foster City
CA) with 50 cycles (94C/60", 60C/3', 72C/6') in a MTC 200

thermo-cycler (MJ Research, Watertown, MA). Each slide was digested twice using 1.3 units of Pvu II (Stratagene) in NEB 2 buffer (New England Biolabs) for 45 minutes at 37 C. After digestion, the tubes were placed at 100 C in 10mM Tris buffer pH 8.0 for 2x5', then blocked with filtered (Millex GV4, Millipore) 1 mg/ml BSA in 2xSSC buffer for 30' at room temperature and rinsed first in 2xSSC 0.1% SDS buffer then in 5xSSC buffer. Each slide was incubated over night at room temperature with a 5xSSC/0.1% Tween 20 buffer containing 1 μ M of the sequencing primer p181 (CGACAGCCGGAAGGAAGAGGGAGC) SEQ ID NO:18 overnight at room temperature. Controls without primer were kept in 5xSSC 0.1% Tween 20 buffer. Glass slides were washed 2 times in 5xSSC 0.1% SDS at 37C for 5' and rinsed in 5xSSC. Primer p181 can hybridize to template B' and the sequence following p181 is CAGCT.... In order to facilitate focusing, green fluorescent beads have been adsorbed to the bottom of the well by incubating each well with 20 μ l of a 1/2000 dilution of 200 nm yellow/green fluorescent, streptavidin coated FluoSpheres^(R) (Molecular Probes, Eugene, OR) in 5X SSC for 20" at room temperature.--

Please replace the paragraph beginning at page 71, line 19, with the following rewritten paragraph:

-- cDNA synthesis - Synthetic mRNA was mixed with mouse liver poly A+ mRNA at different molar ratios (1:1, 1:10, 1:100) and cDNA synthesis on the mixture of synthetic and mouse liver mRNA was performed using the "SMART PCR cDNA synthesis kit" (Clontech, Palo Alto CA) with some minor

modifications. In a cDNA reaction, approximately 1 µg of the mRNA mixture was mixed with the -primer CP5, having at the 5' -end the sequence of CP,

(5'p-AGAAGGAGAAGGAAAGGGAAAGGGTTTTTTTTTTTTTTNN) SEQ ID NO:9.

This primer has been used to make the 1st strand cDNA synthesis. For the 2nd strand synthesis, the "SMART" technique has been used. The basis of the SMART synthesis is the property of the Moloney murine viral reverse transcriptase to add three to five deoxycytosine residues at the 3'-termini of first strand cDNA, when the mRNA contains a 5'-methylguanosine-cap (SMART user manual, Clontech, Palo Alto CA). A CP6 primer, which contains the sequence of CP plus AAAGGGGG (SEQ ID NO:21) at the 3' end,

(5'p-AGAAGGAGAAGGAAAGGGAAAGGGGG) SEQ ID NO:10 has been used for the 2nd strand cDNA synthesis. Buffer and SUPERSCRIPT™ II RNase H- reverse transcriptase from Moloney murine leukemia virus (Life Technologies, Ltd.) were used as described in the instructions and the reaction was carried out at 42 C for 1 hr. The cDNA was assayed by PCR using the primer p251, which contains a fragment of the CP sequence,

(5'-GAGAAGGAAAGGGAAAGG) SEQ ID NO:19 with Taq DNA polymerase (Platinum Taq, Life Technologies, Ltd.).-

Please replace the paragraph beginning at page 72, line 10, with the following rewritten paragraph:

-- Preparation of DNA colonies - The 5'p-cDNA was mixed with different concentrations of the solid phase colony primer, CP2 (5'p-TTTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID

NO:6 and chemically bound to Nucleolink PCR tubes (NUNC) following manufacturer instructions. DNA colonies were then generated using Taq polymerase (AmpliTaq Gold, PE-Applied Biosystems Inc., Foster City CA) with 30 cycles (94C/30", 65C/1', 72C/ 1.5') in a MTC 200 thermo-cycler (MJ Research, Watertown, MA).

Please replace the paragraph beginning at page 74, line 16, with the following rewritten paragraph:

-- Oligonucleotide primers were attached onto Nucleolink plastic microtitre wells (Nunc, Denmark) in order to determine optimal coupling times and chemistries. Oligonucleotides; CP2 (5'-(phosphate)-TTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6, CP8 (5'-(amino-hexamethylene)-TTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6, CP9 (5'(hydroxyl)-TTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6, CP10 (5'-(dimethoxytrityl)-TTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6 and CP11 (5'(biotin)-TTTTTTTTTTAGAAAGGAGAAGGAAAGGGAAAGGG) SEQ ID NO:6, (Microsynth GmbH, Switzerland), were attached to Nucleolink microtitre wells as follows (8 wells each); to each well 20 µl of a solution containing 0.1 µM oligonucleotide, 10mM 1-methyl-imidazole (pH 7.0) (Sigma Chemicals) and 10mM 1-ethyl-3-(3-dimethylaminopropyl)-carbodiimide (pH 7.0) (Sigma Chemicals) in 10mM 1-methyl-imidazole. The wells were then sealed and incubated 50°C for varying amounts of time. The coupling reaction was terminated at specific times by rinsing twice with 200 µl of RS (0.4 N NaOH, 0.25% Tween 20) and twice

with 200 µl TNT (100mM TrisHCl pH 7.5, 150 mM NaCl, 0.1% Tween 20). Tubes were dried at 50°C for 30' and were stored in a sealed plastic bag at 4°C.--

Please replace the paragraph beginning at page 75, line 6, with the following rewritten paragraph:

--Stability was tested under colony growing conditions by adding a PCR mix (20 µl of four dNTPs (0.2 mM), 0.1% BSA, 0.1% Tween 20, 8% DMSO (dimethylsulfoxide, Fluka, Switzerland), IX PCR buffer). The wells were then placed in the thermocycler and for 33 repetitions under the following conditions: 94°C. for 45 seconds, 60°C for 4 minutes, 72°C for 4 minutes. After completion of this program, the wells were rinsed with 5xSSC, 0.1% Tween 20 and kept at 8°C until further use. Prior to hybridization wells are filled with 50 µl 5xSSC, 0.1% Tween 20 heated at 94°C for 5 minutes and stored at RT.

Probe: Oligonucleotide probes, R57 (5'(phosphate)- GTTTGGGTTGGTTGGGTTGGTG (SEQ ID NO:22), control probe) and R58 (5'-(phosphate)-CCCTTCCCTTTCCTTCTCCTTCT (SEQ ID NO:23), which is complementary to CP2, CP8, CP9, CP10 and CP11) were enzymatically labeled at their 5' end terminus with [³²P]dATP (Amersham, UK) using the bacteriophage T4 polynucleotide kinase (New England Biolabs, Beverly, MA). Excess ³²P dATP was removed with a Chroma Spin column TE-10 (Clontech, Palo Alto CA). Radiolabeled oligonucleotides (0.5 µM in 5xSSC, 0.1% Tween 20) were then hybridized to the oligonucleotide derivatized Nucleolink wells at 37°C for two hours. The wells were washed 4 times with 5xSSC, 0.1% Tween 20 at room

temperature, followed by a wash with 0.5xSSC, 0.1% Tween 20 for 15' at 37°C. Wells were then assayed for bound probe by scintillation counting.--

IN THE SEQUENCE LISTING

Please substitute the paper copy Sequence Listing attached hereto for the Sequence Listing originally filed.